

### 68073 Head Rotation Platform for Rats and Mice

The height of the platform can be adjusted according to the specific experimental requirements; the head holder can be adjusted in the vertical direction and 360° rotation in the horizontal direction, the platform below the head is suspended, especially for rats and mice that require 180 °rotation in the head and body during surgery, such as, the experimental study on cerebral ischemia by thread embolism, real-time monitoring the changes of the blood flow in animals' head(fixed fiber on the skull surface) to determine whether the model is successful or not.

68071 Animal rotation platform - Mouse (with mouse rotation adaptor and 60° ear bar)

68074 Animal rotation platform - Rat (with rat rotation adaptor and 18° ear bar)

68073 Animal rotation platform - Rat/Mouse (with rat and mouse rotating adaptor and ear bar)

